▲ DIETZ GROUP

We define mobility



Full catalogue 2024

Quality allows no room for compromise

When it comes to the quality of our products, we - the members of the DIETZ GROUP team – set ourselves the highest standards. Before any rehabilitation aid leaves our factory, it has to go through nume ous testing stations and usage simulations. Only products that truly meet all specifications of our technicians and engineers are deemed to meet our stringent quality standards. As well as going through our internal testing stations, our products are also put through their pa es by independent institutes. Two groups that definitely benefit from our ambition are our retail partners and the people who use our products because they know they can always count on DIETZ quality products to improve their everyday lives.



ISO 13485

As a manufacturer of medical devices, we are certified in accordance with DIN EN ISO 13485:2016. We are meticulous about implementing all requirements tha we - as a medical device manufacturer - need to meet in relation to quality management and quality assurance within the context of our processes.



This confirms that your product has undergone reprocessing in accordance with the MDR and that a professional hygiene service has been provided in line with the recommendations of the Robert Koch Institute.

All DIETZ GmbH products are developed and produced in accordance with the latest provisions of EC Directives 93/42/EEC & 2007/47/EC or 2006/42/EC.

DIETZ GmbH

Reutäckerstraße 12 76307 Karlsbad, Germany info@dietz-group.de www.dietz-group.de

Printing errors, mistakes, price and product changes reserved. Reproduction, in part or in whole, only with written approval. © DIETZ GmbH, 76307 Karlsbad. Dated 08/2024 V07



Hygiene made easy -We'll take care of it for you!

During reprocessing in line with the MDR, our hygiene centre offers a professional hygiene service that complies with the recommendations











Contents

Adaptive wheelchairs	6	
Lightweight manual wheelchairs	8	
Standard manual wheelchairs	10	
DIETZ Power electric wheelchairs	12	
TAiMA rollators	16	
FAKTO+ rollators	18	
Walking aids	19	
Bathroom & toilet products	20	
XL and XXL products	27	
Options	34	
General terms and conditions	39	

PLEASE NOTE

The product portfolio shown is not equally available in all markets – for further details and country-specific item numbers, please speak to your responsible contact person.

AS[01] Adaptive wheelchair



Figure showing
Vicair Active 02 cushin

Options: see order form	AS	G[01]
	up to 125 kg	up to 150 kg
Colour		ments pearl white / petrol / alt blue / black / metallic garnet red
Seat depth (mm)	360 - 450 in steps of 30	480 - 570 in steps of 30
Seat width (mm)	360 - 540 ir	n steps of 20
Seat height, front (mm)	412.5	- 525
Seat height, rear (mm)	387.5 - 525	
Backrest height (mm)	290 - 500 ir	n steps of 30
Backrest angle (degrees)	Adjustable to -5° / 0° / 5° / 10° / 15°	
Push handle height (mm)	860 – 980, depending on backrest height	(optional additional adjustability of +100)
Total width (mm)	Seat width +200	
Total width, folded (mm)	280 with wheels / 300 with wheels and drum brake	
Total height (mm)	Max. 980 with a bad	ckrest height of 500
Total length with leg rests* (mm)	Seat depth +450	Seat depth +450
Total length without leg rests (mm)	Seat depth +330	Seat depth +330
Min. turning radius (mm)	With a seat depth of 450: 1010	With a seat depth of 570: 1150
Armrest height (mm)	220 -	- 260
Armrest width (mm)	6	0
Armrest length (mm)	250 Desk / 35	O Long usable
Lower leg length, in steps of 15 mm	370 - 510	370 - 510
Rear wheels (diameter in inches)	22 / 2	4 / 26
Front wheels (diameter in cm // inches)	12.5/15/17.5 // 5/6/7 12.5/15/17.5 // 5/6/7	
Tipping stability (degrees)	1()°
Max. load capacity (kg)	125	150
Weight (kg)	From 12.8	From 12.8

^{*} Depending on the rear wheel position

AS[01]RF Adaptive wheelchair



Figure showing
Vicair Active 02 cushion + ADI hackre

Options: see order form		1]RF		
	up to 125 kg	up to 150 kg		
Colour	Frame matt black, side eler atomic umbra / titanium grey / coba	nents pearl white / petrol / alt blue / black / metallic garnet red		
Seat depth (mm)	360 - 450 in steps of 30	480 - 510 in steps of 30		
Seat width (mm)	360* - 540 in steps of 20			
Seat height, front (mm)	412.5	- 525		
Seat height, rear (mm)	387.5	- 525		
Backrest height (mm)	290 - 500 ir	n steps of 30		
Backrest angle (degrees)	Adjustable to -5° ,	/0°/5°/10°/15°		
Push handle height (mm)	860 – 980, depending on backrest height (optional additional adjustability of +100)			
Total width (mm)	Seat width +200			
Total width, folded (mm)	280 with wheels / 300 wi	th wheels and drum brake		
Total height (mm)	Max. 980 with a bac	krest height of 500		
Total length with leg rests** (mm)	With leg rests at 77°: seat depth +550 /	/ with leg rests at 82°: seat depth +520		
Min. turning radius (mm)	With a seat depth of 450: 960	With a seat depth of 510: 1060		
Armrest height (mm)	220 -	- 260		
Armrest width (mm)	6	0		
Armrest length (mm)	250 Desk / 35	0 Long usable		
Lower leg length, in steps of 15 mm	370 - 510	370 - 510		
Rear wheels (diameter in inches)	22/2	4 / 26		
Front wheels (diameter in cm // inches)	12.5/15 // 5/6 12.5/15 // 5/6			
Tipping stability (degrees)	10°			
Max. load capacity (kg)	125	150		
Weight (kg)	From 12	From 12.5		

^{*}Seat width of 360 only available without cranked frame ** Depending on rear wheel position

TOMTAR MR-LG

Lightweight folding wheelchair



Options on p. 34			TOMTAR MR-LG		
Colour	Metallic manganese grey				
Seat width (mm)	390	420	450	480	510
Seat depth (mm)	400 / 450	400 / 450	400 / 450	400 / 450	400 / 450
Seat height (mm)	460 / 485 / 510	460 / 485 / 510	460 / 485 / 510	460 / 485 / 510	460 / 485 / 510
Backrest height (mm), adjustable	390 - 480	390 - 480	390 - 480	390 - 480	390 - 480
Side element height (mm)	220 / 235 / 250	220 / 235 / 250	220 / 235 / 250	220 / 235 / 250	220 / 235 / 250
Total dimensions (W x H x L) (mm)	590 x 945 - 1030 x 1050 - 1110	620 x 945 - 1030 x 1050 - 1110	650 x 945 - 1030 x 1050 - 1110	680 x 945 - 1030 x 1050 - 1110	710 x 945 – 1030 x 1050 – 1110
Total length without leg rests (mm)	800	800	800	800	800
Size when folded (mm)	310	310	310	310	310
Max. load capacity (kg)	130	130	130	130	130
Weight* (kg)	15.2	15.3	15.8	16.3	16.5
DE item no.	109039	109042	109045	109048	109051
DE item no. Drum brake	109139	109142	109145	109148	109151

^{*} Complete with leg rests and side elements, without drum brake

CANEO_S Lightweight folding wheelchair



Options on p. 34			CAN	EO_S		
Colour	Ivory	Ivory	Ivory	Ivory	Ivory	Ivory
Seat width (mm)	360	390	420	450	480	510
Seat depth (mm)	400 / 440	400 / 440	400 / 440	400 / 440	400/440	400 / 440
Seat height (mm)	420 - 510	420 - 510	420 - 510	420 - 510	420 - 510	420 - 510
Backrest height (mm)	400 / 430	400 / 430	400 / 430	400 / 430	400 / 430	400 / 430
Side element height (mm)	240 / 260 / 280	240/260/280	240/260/280	240/260/280	240/260/280	240/260/280
Total dimensions (W* x H x L) (mm)	560 x 830 - 920 x 1070 - 1120	590 x 830 - 920 x 1070 - 1120	620 x 830 - 920 x 1070 - 1120	650 x 830 - 920 x 1070 - 1120	680 x 830 - 920 x 1070 - 1120	710 x 830 - 920 x 1070 - 1120
Total length without leg rests (mm)	780	780	780	780	780	780
Size when folded (mm)	310	310	310	310	310	310
Max. load capacity*** (kg)	125	125	125	125	130	130
Weight** (kg)	16.3	16.6	16.9	17.3	18	18.2
DE item no.	210136	210139	210142	210145	210148	210151
DE item no. Drum brake	211136	211139	211142	211145	211148	211151

 $^{^{\}star}$ Plus 10 mm with drum brake ** Complete with leg rests and side elements *** With double cross brace

CANEO_L Lightweight folding wheelchair



Options on p. 34		CANEO_L	
Colour	Ivory	Ivory	Ivory
Seat width (mm)	450	480	510
Seat depth (mm)	470 / 510	470 / 510	470 / 510
Seat height (mm)	480 - 570	480 - 570	480 - 570
Backrest height (mm)	460	460	460
Side element height (mm)	240 / 260 / 280	240 / 260 / 280	240 / 260 / 280
Total dimensions (W* x H x L) (mm)	650 x 910 - 1005 x 1130 - 1200	650 x 910 - 1005 x 1130 - 1200	650 x 910 - 1005 x 1130 - 1200
Total length without leg rests (mm)	850	850	850
Size when folded (mm)	310	310	310
Max. load capacity (kg)	150	150	150
Weight** (kg)	18.7	19	19.3
DE item no.	222145	222148	222151
DE item no. Drum brake	223145	223148	223151

 $^{^\}star$ Plus 10 mm with drum brake $\,^{\star\star}$ Complete with leg rests and side elements

TOMTAR MR-S Standard folding wheelchair



Options on p. 34			TOMTAR MR-S		
Colour	Metallic manganese grey				
Seat width (mm)	390	420	450	480	510
Seat depth (mm)	400 / 450	400 / 450	400 / 450	400 / 450	400 / 450
Seat height (mm)	460 / 485 / 510	460 / 485 / 510	460 / 485 / 510	460 / 485 / 510	460 / 485 / 510
Backrest height (mm)	430	430	430	430	430
Side element height (mm)	220 / 235 / 250	220 / 235 / 250	220 / 235 / 250	220 / 235 / 250	220 / 235 / 250
Total dimensions (W* x H x L) (mm)	590 x 945 x 1050 - 1110	620 x 945 x 1050 - 1110	650 x 945 x 1050 - 1110	680 x 945 x 1050 - 1110	710 x 945 x 1050 - 1110
Total length without leg rests (mm)	800	800	800	800	800
Size when folded (mm)	310	310	310	310	310
Max. load capacity (kg)	130	130	130	130	130
Weight* (kg)	18.1	18.3	18.5	19.3	20.3
DE item no.	108039	108042	108045	108048	108051
DE item no. Drum brake	108139	108142	108145	108148	108151

^{*} Complete with leg rests and side elements, without drum brake

Standard manual wheelchairs


CANEO_B
Standard folding wheelchair



Options on p. 34		
Colour	Anthracite	Anthracite
Seat width (mm)	450	480
Seat depth (mm)	400 / 440	400 / 440
Seat height (mm)	480 / 510	480 / 510
Backrest height (mm)	400	400
Side element height (mm)	230	230
Total dimensions (W x H x L) (mm)	650 x 890 - 920 x 1070 - 1120	680 x 890 - 920 x 1070 - 1120
Total length without leg rests (mm)	780	780
Size when folded (mm)	310	310
Max. load capacity** (kg)	130	140
Weight* (kg)	18.8	19.5
DE item no.	110045	110048

 $^{^\}star$ Complete with leg rests and side elements $\,\,^{\star\star}$ With double cross brace

DECprotect Air cell cushions





DECprotect

Number of chambers	1	1	2	2
Height (mm)	100	60	100	60
Max. load capacity (kg)	150	150	150	150
Weight (kg)	Approx. 0.9	Approx. 0.9	Approx. 0.9	Approx. 0.9
Dimensions W x D (mm)	360 x 360	360 x 360	360 x 360	360 x 360
DE item no.	420027	420033	420020	420039
Dimensions W x D (mm)	360 x 400	360 x 400	360 x 400	
DE item no.	420026	420032	420019	
Dimensions W x D (mm)	400 x 400	400 x 400	400 x 400	400 x 400
DE item no.	420028	420034	420021	420040
Dimensions W x D (mm)	440 x 400	440 x 400	440 x 400	440 x 400
DE item no.	420029	420035	420022	420041
Dimensions W x D (mm)	460 x 460	460 x 460		
DE item no.	420068	420067		
Dimensions W x D (mm)	480 x 440	480 x 440	480 x 440	480 x 440
DE item no.	420030	420036	420023	420042
Dimensions W x D (mm)	480 x 480		480 x 480	
DE item no.	420031		420024	



DIETZ Power electric wheelchairs

SANGO advanced SEGO comfort seat system Electric wheelchair







		RWD	FWD	MWD
The specifications apply to standard I	FWD/RWD/MWD configuratio	ns – for further configura	itions, please refer to the order f	orm.
Colour		Matt gold / I	matt blue / metallic red / gloss	black / gloss white
Max. speed (km/h)			6 / 10 / 12.5	
Max. user weight (kg)	160			
Total weight incl. batteries (kg)	From 154.5			
Batteries (Ah)	60 GEL / 76 GEL / 100 AGM Super Cycle			
Total length with leg rests (mm)	960 - 1760			
Total width* (mm)	615 - 700			
Total height without head rest* (mm)			1010 - 1170	
Seat height, front edge* (mm)			390 - 480	
Effective seat width (mm)			420 - 540	
Effective seat depth (mm)			390 - 560	
Seat tilt (degrees)	Manual // Electr.*		0° / 4° / 8° // -5° - 45° or -5°	- 52°
Backrest angle (degrees)	Manual // Electr.*		90° - 120° // 90° - 120°)
DE item no.		637037	637036	637038

SANGO slimline SEGO comfort seat system Electric wheelchair



		RWD	FWD	MWD	
The specifications apply to standard F	WD/RWD/MWD configuration	ons – for further configu	urations, please refer to the order for	m.	
Colour		Matt gold / matt blue / metallic red / gloss black / gloss white			
Max. speed (km/h)			6 / 10 / 12.5		
Max. user weight (kg)			136		
Total weight incl. batteries (kg)			From 135		
Batteries (Ah)		60 AGM Super Cycle, compact			
Total length with leg rests (mm)		755 - 1760			
Total width* (mm)			530 - 680		
Total height without head rest* (mm)			1030 - 1190		
Seat height, front edge* (mm)			410 - 500		
Effective seat width (mm)			420 - 540		
Seat depth (mm)			420 - 560		
Seat tilt (degrees)	Manual // Electr.*		0° / 4° / 8° // -5° - 45° or -5° -	52°	
Backrest angle (degrees)	Manual // Electr.*		90° - 120° // 90° - 120°		
DE item no.		637046	637045	637047	

^{*}Depending on the model, electrical options such as leg rests with adjustable angles and lengths, an uplifting seat $(0-25^\circ)$ and a seat tilt mechanism (45°) 52 $^\circ$) can be added.

SANGO advanced / lift SEGO comfort seat system Electric wheelchair



		RWD	FWD	MWD	
The specifications apply to standard	FWD/RWD/MWD configurati	ions – for further configu	rations, please refer to the order form		
Colour		Matt gold ,	/ matt blue / metallic red / gloss blac	k / gloss white	
Max. speed (km/h)		6 / 10 / 12.5			
Max. user weight (kg)			160		
Total weight incl. batteries (kg)		From 181			
Batteries (Ah)		60 GEL / 76 GEL / 100 AGM Super Cycle			
Total length with leg rests (mm)		960 - 1760			
Total width* (mm)		615 - 700			
Total height without head rest* (mm)			1010 - 1170		
Seat height, front edge* (mm)		390 - 480			
Electric seat lift (mm)			Up to 300		
Effective seat width (mm)			420 - 540		
Effective seat depth (mm)			390 - 560		
Seat tilt (degrees)	Manual // Electr.*		0° / 4° / 8° // -5° - 45° or -5° - 52	00	
Backrest angle (degrees)	Manual // Electr.*		90° - 120° // 90° - 120°		
DE item no.		637040	637039	637041	

SANGO slimline / lift SEGO comfort seat system Electric wheelchair



		RWD	FWD	MWD	
The specifications apply to standard	FWD/RWD/MWD configurat	ions – for further confi	gurations, please refer to the order form.		
Colour		Matt gol	d / matt blue / metallic red / gloss blac	k / gloss white	
Max. speed (km/h)			6 / 10 / 12.5		
Max. user weight (kg)			136		
Total weight incl. batteries (kg)			From 161.5		
Batteries (Ah)			60 AGM Super Cycle, compact		
Total length with leg rests (mm)		755 - 1760			
Total width* (mm)		530 - 680			
Total height without head rest* (mm)		1030 - 1190			
Seat height, front edge* (mm)		410 – 500			
Electric seat lift (mm)		Up to 300			
Effective seat width (mm)		420 - 540			
Effective seat depth (mm)			420 - 560		
Seat tilt (degrees)	Manual // Electr.*	0° / 4° / 8° // -5° - 45° or -5° - 52°			
Backrest angle (degrees)	Manual // Electr.*	90° – 120° // 90° – 120°			
DE item no.		637049	637048	637050	

SANGO slimline SEGO junior seat system Electric wheelchair



		RWD	FWD	MWD		
The specifications apply to the standard of	specifications apply to the standard configurations – for further configurations, please refer to the order form.					
Colour		Matt gold / n	natt blue / metallic red / gloss bla	ck / gloss white		
Max. speed (km/h)			6 / 10 / 12.5			
Max. user weight (kg)			75			
Total weight incl. batteries (kg)			From 125			
Batteries (Ah)			60 AGM Super Cycle, compact			
Total length with leg rests (mm)		980 - 1410				
Total width* (mm)		530 - 700				
Total height without head rest* (mm)		930 - 1320				
Seat height, front edge* (mm)		410 - 480				
Electric seat lift (mm)		Up to 300				
Effective seat width (mm)		320 - 370				
Seat depth (mm)		300 - 400**				
Seat tilt (degrees)	Manual // Electr.*	0° / 4° / 8° // 0° - 45° or 0° - 52°				
Backrest angle (degrees)	Manual // Electr.*	90° - 120° // 90° - 120°				
DE item no.		637052 637051 637053				

^{*} Depending on the model, electrical options such as leg rests with adjustable angles and lengths, a lift and a seat tilt mechanism (45°/52°) can be added.

** Additional seat depth of 90 mm in combination with lift and seat tilt.

MINKO Electric wheelchair





MINKO

DE item no. (for designs with elect	rical options)	377300 377301 377302		
DE item no.		630000 630001 630002		
Backrest angle (degrees)	Manual // Electr.*	-10° / 0° / 10° / 20° // 90° - 120°		
Seat tilt (degrees)	Manual // Electr.*	-7° to +14°** // 0° - 25°		
Seat depth (mm)		410 - 515		
Nominal seat width (mm)		360 410 460		
Seat height* (mm)		425 - 475		
Total height without head rest* (n	nm)	970 - 1020		
Total width* (mm)		610 - 640		
Total length with leg rests (mm)		1070		
Batteries (Ah)			50 AGM / 60 GEL	
Total weight incl. 45 Ah batteries	(kg)		96.5	
Max. user weight (kg)			136	
Max. speed (km/h)		6		
Colour		Anthracite		
The specifications apply to the star	ndard configurations – for further c	onfigurations, please refe	r to the order form.	

^{*}Depending on the model, electrical options such as leg rests with adjustable angles and lengths, an uplifting seat (0 - 25°), a lift and a seat tilt mechanism (45°/52°) can be added.
**Dependent on the seat height.



TAiMA Ventum

Lightweight rollator







Options on p. 36	TAiMA Ventum		
Colour	Metallic amarena	Metallic cashmere	
Height (mm)	775 – 950	775 - 950	
Seat height x seat width (mm)	610 x 420	610 x 420	
Seat depth (mm)	200	200	
Length x width (mm)	680 x 625	680 x 625	
Width when folded (mm)	240	240	
Max. load capacity (kg)	150	150	
Handle height, adjustable (mm)	7 x 25	7 x 25	
Wheel size (mm)	220 x 32	220 x 32	
Max. additional load in the form of a bag (kg)	5	5	
Weight (kg)	7.7* / 8.2**	7.7* / 8.2**	
DE item no.	200165	200166	

 $^{^{\}star}$ Without / ** With bag, back support and cane holder

TAIMA M-GT

Lightweight rollator, frame size M

TAIMA S-GT

Lightweight rollator, frame size S



Options on p. 36	TAIMA M-GT	TAIMA S-GT
Colour	Metallic reseda green	Metallic sepia
Height (mm)	790 – 950	760 - 860
Seat height x seat width (mm)	590 x 450	510 x 390
Length x width (mm)	680 x 630	620 x 580
Width when folded (mm)	230	230
Max. load capacity (kg)	150	150
Handle height, adjustable	5 x 40 mm	5 x 25 mm
Wheel size (mm)	200 x 35	200 x 35
Max. additional load in the form of a bag (kg)	5	5
Weight (kg)	5.5* / 5.8**	5.3* / 5.6**
DE item no.	200160	200153

 $^{^{\}star}$ Without / ** With bag, back support and cane holder

TAIMA XC

Lightweight rollator with pneumatic tyres

TAIMA XC PU

Lightweight rollator with PU tyres





Options on p. 36	TAIMA XC	TAIMA XC PU
Colour	Metallic red	Metallic oasis green
Height (mm)	850 - 950	850 - 950
Seat height x seat width (mm)	610 x 450	610 x 450
Length x width (mm)	750 x 670	750 x 650
Width when folded (mm)	270	250
Max. load capacity (kg)	150	150
Handle height, adjustable	5 x 25 mm	5 x 25 mm
Wheel size (mm)	310 x 50	310 x 50
Max. additional load in form of a bag (kg)	5	5
Weight (kg)	9.1	9.8
DE item no.	200107	200106

TAIMA M Lightweight rollator

TAIMA M-ECO

Lightweight rollator





Options on p. 36	TAIMA M	TAIMA M-ECO
Colour	Metallic brown	Metallic anthracite
Height (mm)	790 - 950	790 - 950
Seat height x seat width (mm)	590 x 450	590 x 450
Length x width (mm)	680 x 630	680 x 630
Width when folded (mm)	230	230
Max. load capacity (kg)	150	150
Handle height, adjustable	5 x 40 mm	5 x 40 mm
Wheel size (mm)	200 x 35	200 x 35
Max. additional load in the form of a bag (kg)	5	5
Weight (kg)	6.5* / 6.8**	7.8* / 8.1**
DE item no.	200059	200061

 $^{^{\}star}$ Without / ** With bag, back support and cane holder

FAKTO+ rollators

FAKTO+ Standard rollator

FAKTO+ Arthritis

Rollator with armrests





Options on p. 36	FAKTO+	FAKTO+ Arthritis
Colour	Metallic manganese grey	Metallic manganese grey
Total length (mm)	640	620
Total width (mm)	580	580
Total height (mm)	810 - 950	1155 - 1325
Height up to armrests (mm)	-	890 - 1145
Dimensions when folded and lying flat (mm)	910 x 580 x 230	1120 x 580 x 505
Seat width (mm)	380	380
Seat depth (mm)	190	190
Seat height (mm)	595	595
Wheel size (mm)	200 x 50	200 x 50
Max. load capacity (kg)	130	120
Weight incl. cane holder (kg)	8.3*	10*
DE item no.	346NPU	347NPU

^{*} Without basket and tray

Walking aids

TOMTAR G-2E Rigid/reciprocal walking frame



	TOMTAR G-2E	
Description	Walking frame, rigid	
Colour	Black / metallic grey	
Min. total dimensions L x W x H (mm)	460 x 650 x 760	
Dimensions when folded and standing upright L x W x H (mm)	150 x 490 x 760	
Height (mm)	760 - 925	
Recommended body height (cm)	155 - 180	
Max. load capacity (kg)	130	
Weight (kg)	2.8	
DE item no.	200315	

TOMTAR VFG

Quad cane



TOMTAR VFG

DE item no.	200380		
Weight (kg)	0.9		
Max. load capacity (kg)	100		
Recommended body height (cm)	130 - 180		
Ground clearance (mm)	140		
Base L x W x H (mm)	340 / 290 / 210		
Total dimensions L x W x H (mm)	340 x 290 x 720 - 920		
Handle	Ergonomic, for left-handed or right-handed use		
Colour	Black / metallic grey		
Description	Quad cane		

KANJO Bath lift





Options on p. 37	KANJO B-protect	KANJO Power/Eco	
Colour	White / aquamarine	White (optional: blue cover)	
Seat width (mm)	375 / 710 (with side flaps)	375 / 710 (with side flaps)	
Seat depth (mm)	500	500	
Seat height (mm)	65 - 455	65 - 455	
Backrest width (mm)	350	350	
Backrest height (mm)	680	680	
Footrest L x W (mm)	530 x 220	530 x 220	
Max. load capacity (kg)	140	140	
Weight incl. manual control (kg)	10.9	10.8	
DE item no.	750/30	Power: 750/50 (white) 750/62 (blue)	
DE item no.		Eco: 750/51 (white) 750/63 (blue)	

TAYO

Shower stool, round










	TAYO B-protect	TAYO B-protect	TAY0	TAY0
Version	Round	Round, with swivel plate and seat pad	Round	Round, with swivel plate and seat pad
Colour	Light grey / black	Light grey / black	White / matt silver	White / matt silver
Seat width (mm)	320	320	320	320
Seat depth (mm)	320	320	320	320
Seat height (mm)	435 - 585	445 - 595	435 - 585	445 - 595
Total dimensions W x L (mm)	480 - 540 x 430 - 480	480 - 540 x 430 - 480	480 - 540 x 430 - 480	480 - 540 x 430 - 480
Max. load capacity (kg)	130	130	130	130
Weight (kg)	2.1	2.7	2.1	2.7
DE item no.	430116	430117	430106	430107

TAYO











Options on p. 37	TAYO B-protect Shower chair	TAYO Shower stool	TAYO Shower stool	TAYO Shower chair
Version	With armrests	Without armrests	With armrests	With armrests
Colour	Light grey / black	White / matt silver	White / matt silver	White / matt silver
Seat width (mm)	420	500	420	420
Seat depth (mm)	420	420	420	420
Seat height (mm)	440-590	440-590	440-590	440-590
Backrest height (mm)	370	-	-	370
Total dimensions W x L x H (mm)	530-545 x 500-520 x 810-960	530-545 x 500-520 x 440-590	530-545 x 500-520 x 620-780	530-545 x 500-520 x 810-960
Max. load capacity (kg)	130	130	130	130
Weight without accessories (kg)	4.58	3.3	3.92	4.58
DE item no.	430114	430101	430102	430104

TAYO

Shower stool, triangular





	TAYO B-protect shower stool	TAYO shower stool
Version	Triangular	Triangular
Colour	Light grey / black	White / matt silver
Seat width (mm)	380	380
Seat depth (mm)	340	340
Seat height (mm)	440 - 590	440 - 590
Total dimensions W x L (mm)	430 - 450 x 395 - 410	430 - 450 x 395 - 410
Max. load capacity (kg)	130	130
Weight (kg)	2.1	2.1
DE item no.	430115	430105





22 DIETZ GROUP // Full catalogue DIETZ GROUP // Full catalogue 23

SmartFix Toilet seat raiser







SmartFix B-protect	SmartFix	SmartFix
With armrests	With armrests	Without armrests
White	White	White
80 / 120 / 150	80 / 120 / 150	80 / 120 / 150
200	200	-
255 x 195	255 x 195	255 x 195
560 x 497 x 377	560 x 497 x 377	455 x 497 x 177
120	120	120
60 each	60 each	-
4.2	4.2	2.82
429114	429112	429113
	With armrests White 80 / 120 / 150 200 255 x 195 560 x 497 x 377 120 60 each 4.2	With armrests With armrests White White 80 / 120 / 150 80 / 120 / 150 200 200 255 x 195 255 x 195 560 x 497 x 377 560 x 497 x 377 120 120 60 each 60 each 4.2 4.2

Swivel aid For bath lift

Swivel and transfer aid

For bath lift





	Swivel aid	Swivel and transfer aid
Colour	White / with blue cover	White / with blue cover
Total width (mm)	320	320
Total length (mm)	320	411
Total height (mm)	22	22
Diameter (mm)	320	320
Max. load capacity (kg)	140	140
Weight (kg)	1.1	1.3
DE item no.	750/12	750/14

Relaxon

Toilet seat raiser







Options on p. 37	Relaxon B-protect	Relaxon	Relaxon
Version	With lid	With lid	Without lid
Colour	White	White	White
Seat height (mm)	100	100	100
Opening L x W (mm)	277 x 205	277 x 205	277 x 205
Total dimensions W x L x H (mm)	370 x 400 x 100	370 x 400 x 100	370 x 400 x 100
Max. load capacity (kg)	200	200	200
Weight (kg)	1.4	1.4	1
DE item no.	430302	000111	000110

Santilo

Bath board with handle



	San	tilo
Description	Bath board with handle	Bath board with handle
Colour	White	White
Adjustment range (mm)	300 - 530	380 - 650
Seat surface W x D (mm)	510 x 240	590 x 240
Total width (mm)	350	350
Total length (mm)	690	740
Total height (mm)	160	160
Max. load capacity (kg)	150	150
Weight (kg)	2.2	2.7
DE item no.	140000	140001





24 DIETZ GROUP // Full catalogue DIETZ GROUP // Full catalogue 25

Levina 320-2

Shower / commode transport wheelchair

Levina 400-2

Shower / commode wheelchair







	Levina 320-2	Levina 400-2
Description	Shower / commode transport wheelchair	Shower / commode wheelchair
Colour	White	White
Tyres	4 x 5"	2 x 5", 2 x 24"
Seat dimensions (W x D x H) (mm)	435 x 425 x 530	435 x 425 x 530
Backrest height (mm)	400	400
Total dimensions W x L x H (mm)	555 x 885 - 925 x 930	635 x 1015 - 1055 x 945
Drive-over height (mm)	435	435
Max. load capacity (kg)	120	120
Weight (kg)	11.5	17
DE item no.	140021	140022

9400-A/9401-A/ 9402-A/9404-A Shower chair, shower stool

and folding shower seat









120 4.1	120 3.9	120 2.75	120
	120	120	
020 % 000 % 000			170 X 100 X 170 070
525 x 530 x 880 - 950	525 x 520 x 860 - 930	535 x 440 x 640 - 710	470 x 490 x 470 - 570
180	180	180	-
460 - 530	460 - 530	460 - 530	470 - 570
400	400	400	400
400	400	400	400
White	White	White	White
Shower chair, folding	Shower chair	Shower stool	Shower seat, folding
9400-A	9401-A	9402-A	9404-A
	Shower chair, folding White 400 400 460 - 530	Shower chair, folding Shower chair White White 400 400 400 400 460 - 530 460 - 530 180 180	Shower chair, folding Shower chair Shower stool White White White 400 400 400 400 400 400 460 - 530 460 - 530 460 - 530 180 180 180

TSU-2/TSU-4 Universal commode wheelchair



Options on p. 37	TSU-2*	TSU-4*
Description	Commode wheelchair	Commode wheelchair
Colour	Metallic manganese grey	Metallic manganese grey
Seat dimensions (W x D x H) (mm)	470 x 440 x 510	470 x 440 x 510
Side element height (mm)	270	270
Drive-over height (mm)	420	420
Total dimensions W x L x H (mm)	560 x 815 - 825 x 980	560 x 815 - 825 x 980
Max. load capacity (kg)	120	120
Weight (kg)	12.9	13.1
Wheels	4 x 5" wheels	4 x 5" wheels
DE item no.	401001	401002

^{*} The TSU-2 features parking brakes on the two rear wheels only. The TSU-4 features parking brakes on all four wheels.

TOMTAR DH/TOMTAR DS Shower stool/chair

Comfort Shower stool







	TOMTAR DH	TOMTAR DS	Comfort shower stool
Description	Shower stool	Shower chair	Shower stool
Colour	White	White	White
Seat width (mm)	380	380	430
Seat depth (mm)	295	295	405
Seat height (mm)	Internal 355 – 505	Internal 355 - 505	460 - 560
	External 385 - 535	External 385 - 535	460 - 560
Side element height (mm)	-	-	180
Total dimensions W x L x H (mm)	490 x 400 x 385 - 535	490 x 475 x 680 - 830	558 x 620 x 510 - 540 x 640 - 740
Max. load capacity (kg)	100	100	160
Weight (kg)	2	2.8	3.6
DE item no.	411300	411301	140012

26 DIETZ GROUP // Full catalogue

Bano Swivel bath seat

9403-A/B Bath seat









	Bano	9403-A	9403-B
Description	Swivel bath seat	Bath seat without backrest	Bath seat with backrest
Colour	White	White	White
Seat surface W x D (mm)	400 x 400	400 x 400	400 x 400
Total width (mm)	710 - 870	735 - 890	735 - 890
Total length (mm)	570	420	560
Total height (mm)	575	285	430
Max. load capacity (kg)	120	120	120
Weight (kg)	6.5	2.7	3.6
DE item no.	000118	411225	411224

XL and XXL products

tauron|rsi Folding manual XXL wheelchair



Options: see order form	tauron rsi	
Colour	Anodised aluminium frame: moon grey	
Selectable seat width, front (mm)	600 / 650 / 700 / 750	
Adjustable seat width, rear (mm)	0 to - 100	
Selectable seat depth (mm)	430 / 360 (short side frame) 480 / 390 (long side frame)	
Adjustable seat height (mm)	390 / 430 / 470 (without seat cushion)	
Adjustable push bar position (degrees)	360° in 45° steps	
Total width (mm)	Seat width at front +190	
Size when folded (mm)	330	
Total height (mm)	1030	
Total length with leg rests (mm)	1040 incl. standard leg rests with foot board (short side frame) 1090 incl. standard leg rests with foot board (long side frame)	
Total length without leg rests (mm)	790 (short side frame) 840 (long side frame)	
Armrest height (mm)	300 (without seat cushion)	
Armrest width (mm)	50	
Armrest length (mm)	250	
Lower leg length (mm)	Adjustable: 390 – 460	
Rear wheels (diameter in inches)	24" (rim, push rim and hub: black anodised aluminium)	
Front wheels (diameter in mm // inches)	150 / 200 // 6" / 8"	
Max. load capacity (kg)	250	
Weight (kg)	28	
DE item no.	713008	

28 DIETZ GROUP // Full catalogue DIETZ GROUP // Full catalogue 29

CANEO_XL
Heavy-duty folding wheelchair



Options on p. 38	CANEO_)	(L 170 kg	CANEO_X	(L 200 kg
Colour	Ivory	Ivory	Ivory	Ivory
Seat width (mm)	560	600	560	600
Seat depth (mm)	400 / 440	400 / 440	400 / 440	400 / 440
Seat height (mm)	440 / 470	440 / 470	440 / 470	440 / 470
Backrest height (mm)	430	430	430	430
Side element height (mm)	230	230	230	230
Total dimensions (W* x H x L) (mm)	770 x 920 - 950 x 1150 - 1200	800 x 920 - 950 x 1150 - 1200	770 x 920 - 950 x 1150 - 1200	800 x 920 - 950 x 1150 - 1200
Total length without leg rests (mm)	880	880	880	880
Size when folded (mm)	340	340	340	340
Max. load capacity (kg)	170	170	200	200
Weight** (kg)	24	25	24	25
DE item no.	151161	151162	151056	151060
DE item no. Drum brake	151163	151164	151156	151160

^{*} Plus 1 cm with drum brake ** Complete with leg rests and side elements

SANGO XL & XXL SEGO XL & XXL seat system Electric wheelchair





		SANGO xl FWD	SANGO xxl FWD
The specifications apply to the sta	ndard configurations – fo	or further configurations, please refer to the order	form.
Colour		Matt blue	Matt blue
Max. speed (km/h)		6	6
Max. user weight (kg)		200	250
Total weight with batteries (kg)		From 159.5	From 163.5
Batteries (Ah)		60 GEL / 76 GEL / 100 AGM Super Cycle	76 GEL / 100 AGM Super Cycle
Total length with leg rests (mm)		985 - 1680	990 - 1835
Total width* (mm)		695 – 700	740 - 780
Total height* without head rest (mm)		1125	1110 - 1380
Seat height,* front edge (mm)		455 - 515	455 - 520
Electric seat lift (mm)		-	Up to 250
Effective seat width (mm)		540 - 650	600 - 640
Effective seat depth (mm)		480 - 550	490 - 650
Seat tilt (degrees)	Manual // Electr.*	0°/ 4° / 8° // 0° - 25°	0°/4°/8°//0°-25°
Backrest angle (degrees)	Manual // Electr.*	90° – 120°	90° - 120° // 90° - 135°
DE item no.		637021	630042

Depending on the model, electrical options such as leg rests with adjustable angles and lengths, an uplifting seat (0 - 25°), a lift (xxl only) and a seat tilt mechanism (25°) can be added

v-max⁺

Electric power assist for assisted pushing and braking





	v-max ⁺	
Colour	Silver-black	
Dimensions H x W x D (mm)	262 x 600 - 1010 x 167	
Speed (km/h)	5.5 forwards, 3 in reverse	
Max. climbing capacity (%)	Up to 12	
Range of up to (km)	11.5	
(Rechargeable) batteries (V/Ah)	2 x 12 / 22	
Continuous rating of motor (W)	419	
Total permissible weight (kg)	300	
Weight incl. (rechargeable) batteries (kg)	27.7	
DE item no.	346063	

SOLO+ black series

Powerful wheelchair add-on drive





SOLO ⁺ black series		
Colour	Matt black / neon green	
Dimensions of wheelchair & drive H x W x D (mm)	103x seat width at front + $25x$ 104 for short side frame / 109 for long side frame (incl. leg rests)	
Size when folded (mm)	390	
Speed (km/h)	Variable: 6 (forwards), 3 (in reverse)	
Max. climbing capacity (%)	20	
Range of up to (km)	15	
(Rechargeable) batteries (V/Ah)	28.8 / 16	
Continuous rating of motor (W)	150	
Tyres (inches)	24 / 13 / 8	
Total permissible weight (kg)	250	
Weight incl. (rechargeable) batteries (kg)	24.2	
DE item no.	345001	

30 DIETZ GROUP // Full catalogue DIETZ GROUP // Full catalogue 31

Rollator XXL Heavy-duty rollator

Walking frame XXL Walking aid XXL





Options on p. 38	Rollator XXL	Walking frame XXL
Colour	Silver	Black
Seat W x D x H (mm)	600 x 340 x 540	-
Total width (mm)	710	710
Total height (mm), adjustable	820 - 925	760 – 1010
Total length / folded (mm)	730 / 300	505
Width between handles (mm)	-	520
Max. load capacity (kg)	200	270
Weight (kg)	12.6* / 11.8**	6
DE item no.	340NPU	317

^{*} With shopping basket ** Without shopping basket

Quad cane XXL Walking aid



	Quad cane XXL	
Colour	Black	
Length adjustment (mm)	710 - 930	
Max. load capacity (kg)	270	
Weight (kg)	1.9	
DE item no.	10106	

BURE XL

Walker XL, height electrically adjustable



	BURE XL	
Width (mm)	880	
Length (mm)	780	
Height (mm)	910 - 1310	
Max. load capacity (kg)	240	
Weight (kg)	29	
Armrest/internal dimension (mm)	440 -620	
Wheel size (inches/mm)	5" / 125	
DE item no.	240401	

KANJO XL Bath lift XL



DE item no.	750/70	
Total weight incl. manual control (kg)	11.5	
Max. load capacity (kg)	170	
Backrest W x H (mm)	350 x 680	
Seat height (mm)	65 - 455	
Seat depth (mm)	500	
Seat width (mm)	375 / 710 (without / with side flap)	
Total depth, backrest down (mm)	900	
Total depth, backrest up (mm)	690	
Total height (mm)	1105	
Colour	White / with blue covers	
Description	Bath lift XL	
	KANJO XL	

tauron dts Shower commode chair XXL







Figure showing rectangular toilet pan

Options on p. 38		tauron dts	
Side element version	Standard	Abduction with standing aid	Swivelling
Colour	Manganese grey / lime green	Manganese grey / lime green	Manganese grey / lime green
Seat width (mm)	560/630/700	560/630/700	560/630/700
Seat width, front* (mm)	-	710 / 780 / 850	-
Seat depth (mm)	450 / 510	450 / 510	450 / 510
Seat height (mm)	450 - 600	450 - 600	450 - 600
Side element height (mm)	200	200	200
Side element height, front* (mm)	-	330	-
Drive-over height (mm)	343 - 493	343 - 493	343 - 493
Total dimensions W x L x H (mm)	750 - 932 x 705 - 747 x 860 - 1060	810 - 960 x 770 - 806 x 860 - 1060	750 - 932 x 705 - 747 x 860 - 1060
Max. load capacity (kg)	250	250	250
Weight (kg)	14	14.5	15
DE item no.	720000 + 720100	720000 + 720101	720000 + 720102

 $^{^{\}star}$ Only differs from the rear dimensions in the case of abduction side element with standing aid.

tauron|dtrs Shower commode wheelchair XXL







Options on p. 38		tauron dtrs			
Side element version	Standard	Abduction with standing aid	Swivelling		
Colour	Manganese grey / lime green	Manganese grey / lime green	Manganese grey / lime green		
Seat width (mm)	560/630/700	560/630/700	560/630/700		
Seat width, front* (mm)	-	710 / 780 / 850	-		
Seat depth (mm)	450 / 510	450 / 510	450 / 510		
Seat height (mm)	550 - 600	550 - 600	550 - 600		
Side element height (mm)	200	200	200		
Side element height, front* (mm)	-	330	-		
Drive-over height (mm)	443 - 493	443 - 493	443 - 493		
Total dimensions W x L x H (mm)	730 - 885 x 930 - 937 x 964 - 1063	853 - 993 x 1003 - 1009 x 964 - 1063	730 - 885 x 930 - 937 x 964 - 1063		
Max. load capacity (kg)	250	250	250		
Weight (kg)	16	16.5	17		
DE item no.	720001 + 720103	720001 + 720104	720001 + 720105		

^{*} Only differs from the rear dimensions in the case of abduction side element with standing aid.



Options for manual wheelchairs

			Caneo_S	Caneo_B	Caneo_L	TOMTAR MR-LG	TOMTAR MR-S	
	Description		ပိ	చ	ర	유	유	DE item no.
	Angle-adjustable amputee board	Left				•	•	196098
181		Right				•	•	196097
	Angle-adjustable amputee board,	Left		•	•			See order form
Fa	mounted	Right	•	•	•			See order form
Not shown	Adjustable backrest with 3 eyes	Seat width of 390 mm						190042
	and waist belt opening,	Seat width of 420 mm						190043
	mounted	Seat width of 450 mm		•				190044
		Seat width of 480 mm		•				190045
		Seat width of 510 mm						190046
7	Anti-tipping support	Left	•	•	•	•	•	196236
*		Right	•	•	•	•	•	196237
	Lap belt	Seat width of 390 - 450 mm				•	•	196238
	1	Seat width of 480 – 510 mm				•		196239
MI	Press lap belt with metal clasp		•	•	•			402114
de	BW leg rest, infinitely	Left	•	•	•	•	•	See order form
6	angle-adjustable with gas pressure, seat width 360 – 440 mm	Right	•	•	•	•	•	See order form
	BW leg rest, infinitely	Left	•	•	•	•		See order form
	angle-adjustable with gas pressure, seat width 450 – 540 mm	Right	•	•	•	•	•	See order form
	BW leg rest, seat width of 450 - 600 mm,	Left			•			198010
	removable and angle-adjustable, mounted	Right			•			198016
	Brake lever extension for parking brake,	Left	•	•	-			198065
	mounted	Right	•	•	•			198064
	Brake lever extension	Pair				•	•	196240
=	CardLock, locking system,	Anthracite		•				111045
	incl. 500 mm long lock chain and a card for locking/unlocking, mounted (on right)	Wall bracket		•				197024
	Single-hand drive, mounted	Left	•	•	•			197039
		Right	•	•	•			197040
	24" PU tyres, black anodised aluminium rim stainless steel spokes, connecting axle car without tools – wheelchair remains fully fo	n be removed						

PLEASE NOTE

The item numbers
listed here only apply when
the accessories are purchased
together with the main
product. If they are
ordered subsequently,
the item numbers
may differ.

the item ner may differ.	Description		Caneo_S	Caneo_B	Caneo_L	TOMTAR MR-LG	TOMTAR MR-S	DE item no.
	One-hand parking brake,	Left		•		'		190033
** 3g/	mounted	Right	•	•	•			190034
	Black folding lock		•	•	•	•	•	196259
(Push rim cover	Black/grooved	•	•	•	•	•	196256
(-)	for 24" wheel, silicone	Black/smooth	•	•	•	•	•	196257
T	Infusion holder for wheelchair	Approx. 1000 mm with 2 hooks	•	•	•	•	•	190053
es _p	Head rest complete with clamp bracket (not suitable for use in conjunction with tie-down system!)	Seat width of 400 – 560 mm	•	•	•	•	•	196252
	Tie-down system, mounted,	Long lap belt	•	•				See order form
The state of the s	crash test certified according to ISO 7176-19 (approved for transporting people in a vehicle)	Short lap belt	•	•				See order form
+ 1	Wheelchair set (anti-tipping support, press lap belt, cane holder set)		•	•	•	•	•	404076
-	Height-adjustable side element, Desk + Long combination	Pair		•				198020
Not shown	Height-adjustable side element	Pair				•	•	196244
	Spoke protectors, transparent for 24" wheel	Withfasteners	•	•	•	•	•	198069
b ,	Cane holder with high-profile cup	With screw bracket for 25 mm tube	•	•	•			196255
Not shown	Cane holder, with screw bracket for 25 mr	n tube				•	•	196260
	Therapy tray	Seat width of 390 - 450 mm	•	•	•	•		470060
		Seat width of 480 - 510 mm	•	•	•	•	•	470061

36 DIETZ GROUP // Full catalogue 37

Options for rollators

	Description		TAIMA Ventum	TAIMA M-GT	TAIMA S-GT	TAIMA XC	TAIMA XC PU	TAIMA M	TAIMA M-ECO	FAKT0+	FAKTO+ Arthritis	DE item no.
3322	Cover	Grey with lead strip			•	•	•	•	•			297002
	30701	orey with tead only			_			_	_			237002
3 T	One-hand brake system	Mounted on the left	•	•	•		•	•	•			290105
	(TAiMA)	Mounted on the right	•	•	•		•	•	•			290104
	One-hand brake system,	Left								•	•	290107
	retrofit kit	Right								•	•	290108
-	Umbrella	Red	•	•	•	•	•	•	•	•	•	601270
77		Macchiato	•	•	•	•	•	•	•	•	•	601271
	Rollator bag M	W x H x D: 380 x 175 x 300 mm, volume 20 l, max. load 15 kg		•				•	•			290101
	Rollator bag S	W x H x D: 325 x 175 x 250 mm, volume 14 l, max. load 12 kg			•	•	•					290102
4	Premium rollator bag	Black/amarena (label)	•									290406
	(Fidlock), with shoulder strap	Black/cashmere (label)	•									290407
3	Back support (TAiMA)	Standard version for TAIMA S-GT / M-GT / M-ECO				•	•	•				290007
	Comfort back support (TAiMA)	Padded		•	•	•	•	•	•			290038
	Back strap (Fakto+)	Black								•	•	290081
TAIMA	Back support (TAiMA Ventum)	Cashmere	•									290401
TAJMA		Black	•									290402
TARMA		Amarena	•									290400
	Push handles, long (TAiMA)	Mounted Total height 990 – 1140 mm		•				•	•			290106
	Push handles, long (TAiMA Ventum)	Mounted Total height 910 – 1085 mm	•									290408
	Tray				•	•	•	•	•	•	•	290137

Options for bathroom & toilet products

	Description	KANJO B-protect	KANJO XL	KANJO Power / Eco	TSU-2/TSU-4	TAY0 B-protect	TAYO	DE item no.
00	Cover caps, set of 2, light grey					•		430129
	Cover caps, white						•	430125
II	Armrests, pair, light grey					•		430126
W	Armrests, pair, white						•	430122
	Aquamarine covers, B-protect , full-length backrest cushion	•						477010
Not shown	Warranty extension (KANJO) from 3 to 5 years	•	•	•				751/16
100 m	2 cm height adapter (set) incl. 2 suction cups	•	•	•				075167
	Backrest, light grey					•		430127
7	Backrest, white						•	430123
	Side flap protectors	•	•	•				075166
	Seat riser with hygiene opening, padded				•			490002
#	Seat pad, blue						•	430120



The item numbers
listed here only apply when
the accessories are purchased
together with the main
product. If they are
ordered subsequently,
the item numbers
may differ.





Options for XL and XXL products

PLEASE NOTE The item numbers listed here only apply when the accessories are purchased together with the main product. If they are ordered subsequently, the item numbers may differ.

	CANEO_XL (107/20	Rollator XXL	tauron dts	tauron dtrs	
Description	CAN	Roll	tau	tau	DE item no.
Not shown Cover Grey with lead strip		•			297002
Angle-adjustable Left	•				198000
amputee board, mounted Right	•				198004
Anti-tipping support Mounted, left	•				198066
Mounted, right	•				198005
Press lap belt with metal clasp	•				402114
BW leg rest, seat width of 450 - 600 mm, Left	•				198010
removable and angle-adjustable, mounted Right	•				198016
One-hand parking brake, Left	•				190033
factory mounted Right	•				190034
One-hand brake system, retrofit kit Left		•			290107
Right		•			290108
Round toilet Round			•		720107
pan: W x H: 300 x 160 mm / dia. 300 mm,				•	720118
rectangular: L x W x H: 410 x 300 x 160 mm Rectangular			•		720108
				•	720119
Head rest complete with clamp bracket (not suitable for use in conjunction with tie-down system!) Seat width of up to 560 mm	•				196252
Tie-down system, mounted, crash test certified Long lap belt according to ISO 7176-19 (approved for transporting people in a vehicle)	•				195017
Wheelchair set (anti-tipping support, Seat width of 540 mm press lap belt, cane holder set)	•				404076
VICAIR XXTRA 02 seat cushion, seat width 550 - 750 mm, seat depth 450/500 mm, seat height 100 mm	•				See order form
Cane holder Black		•			500007
Tray Transparent		•			500594
Therapy tray Seat width of 560 mm	•				195045
Seat width of 600 mm	•				195047

General terms and conditions version: 01/2019

§ 1 General information

- 1. Our terms and conditions of business, delivery and payment apply to the exclusion of all others. All side agreements in addition to these terms and conditions, even if accepted by representatives or employees must be confirmed by us in writing in order to be valid. General terms and conditions of the customer only apply if we have expressly agreed to them in writing.
- 2 If any individual terms and conditions are invalid this shall not affect the validity of the remaining terms and conditions.

§ 2 Quotation and conclusion of contract

- 1. Our quotations are non-binding. Orders are only binding on us if we confirm them or fulfil them by sending the goods to the customer. Verbal side agreements shall only be binding if we confirm them in writing.
- 2.Confirmationoftheorderobligatesthecustomertotake delivery of the goods detailed in the order confirmation.
- 3. As a basic principle, orders must be placed in writing (which includes by fax) If an order is placed only verbally transcription errors and any misunderstandings shall be the responsibility of the customer. If a written order confirmation exists, the scope and content of the order shall be taken from there.

§ 3 Delivery and delivery dates

- 1. The goods shall be delivered in the specified versions and packing units. We reserve the right to make partial deliveries.
- 2. If the customer is in arrears with a payment owed to § 6 Data storage us, our obligation to deliver shall be suspended.
- 3. If an agreed delivery date is culpably exceeded, we shall only be deemed in default of delivery once a reasonable grace period has expired, which shall commence as soon as we receive the grace period notification. If we default in delivery for reasons for which we are responsible, the customer shall be entitled to claim - for each complete week of default - compensation of 0.5% of the invoice value of the goods and services affected by the default, up to a maximum total of 5%. The above provision in no way changes the onus of proof to our detriment. Any claims over and above this shall be excluded unless gross negligence, as a minimum, is what led to the default.

§ 4 Force majeure

- 1. In the case of force majeure, examples of which include strikes, lockouts or other events that we are unable to avert despite exercising a level of diligence that is reasonable given the circumstances, the contractual obligations of the parties shall - regardless of whether the event occurs within our own organisation or within a supplier's operations - be suspended for the duration of the disruption and in terms of the extent.
- 2. If any resulting delays exceed a period of six weeks, both contracting parties shall be entitled to withdraw from the contract as far as the affected scope of supply and services is concerned. No other entitlements exist.

§ 5 Prices and terms of conditions of payment

1. The prices specified (in EUR) in our sales documents (catalogues, ordering price lists, etc.) apply at the time of publication in each case. Promotional prices are only valid for the time period specified in each case. We reserve the right to change our prices following the pub-

- been concluded, the agreed price cannot be changed. 2. The prices apply ex works from Germany excluding packaging.
- 3. Our prices are exclusive of statutory VAT. This shall be shown separately on the invoice at the statutory rate that applies on the invoicing date.
- 4 Unless otherwise stated in the order confirmation our invoices are payable within 10 days without any deductions. The time of performance shall be stated on the invoice. The customer shall be deemed in default of payment without the need for a reminder.
- 5. If the payment period is exceeded, we shall be entitled to charge default interest of at least 4% p.a. starting from when the customer defaults. Both we and the customer reserve the right to demonstrate a lower \$ 10 Reservation of title or higher level of loss or damage in individual cases.
- 6. In the event of default in payment and justified doubts concerning the solvency or creditworthiness of the customer, we shall be entitled to demand security or advance payment for any deliveries not yet made. demand the immediate settlement of all claims arising under the business relationship and take back goods delivered under reservation of title.
- 7. The customer shall only be entitled to set off claims if their counterclaims have been established as final and absolute, are undisputed or have been recognised by us.

Notice in accordance with Section 33 of the German Federal Data Protection Act: we store customer data.

§ 7 Shipping and damage in transit

- 1. The goods shall be shipped at the customer's risk. The risk shall pass to the customer as soon as the goods are dispatched from the warehouse. Insurance shall only be provided if expressly requested. In this case, the insurance premium charged shall be equal to 1% of the value of the goods.
- 2. If the customer identifies damage to the packaging (damage in transit), they must get the transport company to make an official note of the damage when taking receipt of the goods. Damage in transit that is only identified once the goods are unpacked must be reported to us in writing and received by us within three days of receipt of the goods.

§ 8 Warranty and liability

- 1. Complaints that relate to defects and to the properties and condition of the delivered goods shall only be considered if they are submitted in writing within eight days of receipt of the goods or, in the case of hidden defects, within 8 days of their discovery - subiect to a maximum of six months after receipt of the goods - along with supporting documents.
- 2. Our statutory warranty obligation shall be limited, at our option, to either rectification of the defects or a substitute delivery. If rectification of the defects/substitute delivery is unsuccessful, the customer shall at their ontion be entitled to rescind the contract or to demand an appropriate reduction in the purchase price.
- 3. We cannot be held liable for consequential loss

- lication of the sales documents. Once a contract has 4. The limitation period for claims and rights arisin from defects is 1 year
 - 5. Limitations and exclusions of liability do not apply in the case of wilful intent or gross negligence.

Goods can only be returned with our consent. The reason for return must be specified on the return shipment request and the relevant copy of the invoice must be enclosed with the return shipment. Special orders cannot be taken back. When goods are returned in exchange for a credit note, a 15% administrative fee and the outward freight costs shall be deducted. All returns must be processed in such a way that no costs are incurred by us. We cannot take back items that have been modified or damaged by the purchaser.

- 1. The purchased goods shall remain our property until we receive full payment of the receivables that are owed to us under the business relationship with the purchaser. The purchaser is permitted to dispose of the nurchased goods in the ordinary course of business The purchaser hereby assigns to us, by way of security. the receivables due from third parties as a result of resale. The purchaser is authorised to collect these for our account until further notice or until their navments to us cease
- 2. In the event that the goods belonging to us are accessed by third parties, particularly bailiffs, the purchaser shall inform them that they are our property and notify us immediately by registered letter or fax.
- 3. In the event of the purchaser acting in breach of the contract, particularly in the case of default in payment, we shall be entitled to take back the goods belonging to us. Exercising our rights under reservation of title does not constitute withdrawal from the contract.

§ 11 Place of performance, place of jurisdiction, applicable law

- 1. The place of performance for all obligations arising under the contractual relationship is Karlsbad, Germany
- 2. As regards all disputes that arise under the contractual relationship, legal action is to be taken at the court with jurisdiction over our registered place of business whenever the customer has "full merchant" status under German law, is a legal person under public law or is a special fund under public law. We are also entitled to take legal action at the location of the customer's head office.
- 3. German law applies exclusively to the exclusion of the UN Convention on Contracts for the International Sale of Goods

§ 12 Copyright

We reserve all title to and copyright in illustrations particularly those showing our products - drawings, costings / calculations and other documents. This especially applies to all those written documents that are identified as "confidential". Before these are used, utilised or passed on, the customer or business partner requires our express written permission

§ 13 Severability clause

If any provision of this contract is ineffective or should become ineffective during the term of this contract, the validity and effectiveness of the other provisions shall remain unaffected

▲ DIETZ GROUP

Adaptive wheelchairs / Manual wheelchairs
Electric wheelchairs by DIETZ Power
Rollators
Walking aids
Bathroom & toilet products
XL, XXL products and custom models

